

# New Mitochondrial Function Assay Technology

Lawrence A. Wiater, Kyle O'Hollaren, Xiang-He Lei, and Barry R. Bochner  
Biolog, Inc. Hayward, CA, USA


## Abstract

We have developed a new mitochondrial function assay technology that measures the rates of metabolism of mitochondrial substrates and the sensitivity of metabolism of these substrates to mitochondrial inhibitors. The technology employs saponin permeabilized cells and a redox dye added to 96-well microplates that contain mitochondrial substrates or inhibitors precoated and dried into the wells. The MitoPlate S-1™ has a triplicate repeat of a set of 31 substrates. Mitochondrial function is assayed by measuring the rates of dye reduction from electrons flowing into and through the electron transport chain from substrates whose oxidation produces NADH (e.g., L-malate) or FADH<sub>2</sub> (e.g., succinate). The electrons donated to complex 1 or complex 2 travel to the distal portion of the electron transport chain where a tetrazolium redox dye (MC) acts as a terminal electron acceptor and changes from colorless to a purple formazan upon reduction. All 96 assays in the MitoPlate are run concurrently, and each assay provides different information because each substrate follows a different metabolic route using different transporters to enter the mitochondria, and then different dehydrogenases to produce NADH or FADH<sub>2</sub>. The MitoPlate S-1™ can also be used to assess the activity and specificity of substrate transport inhibitors, dehydrogenase inhibitors, or electron transport chain inhibitors. A second assay plate, the MitoPlate I-1™, provides another assessment of mitochondrial function by measuring the sensitivity of mitochondrial electron flow to a set of 22 diverse inhibitors titrated at 4 dilutions. The I-1 plates can be run using any of the NADH or FADH<sub>2</sub> producing substrates, each providing additional information. Using these new assays we show that the mitochondria from different cell types exhibit different functional properties. This new technology will assist efforts to understand how mitochondria change in cell models of human disorders that have a mitochondrial basis.

## The Assay Technology


**Figure 1** outlines the simple assay protocol. 30 µl of a permeabilizing assay mix containing saponin and Redox Dye MC in isotonic buffer is pipetted into all wells and incubated at 37° C for 1 hour. To start the assay, 30 ul of a cell suspension in isotonic buffer is added to each well. The recommended cell density is 1,000,000 cells/ml resulting in 30,000 cells/well. To record the rates of dye reduction in the wells, the MitoPlate is loaded into the OmniLog, which reads at 5 minute intervals for 2 to 4 hours. For MitoPlate I-1 with 22 mitochondrial inhibitors, the permeabilizing assay mix also contains an NADH or FADH<sub>2</sub> producing substrate such as L-malate or succinate. **Figure 3 and 4** show, respectively, the test layout in the MitoPlate S-1 and the MitoPlate I-1.

**Figure 1.** The Assay Protocol.



## Assays and Results

**Figure 2.** With this assay technology, mitochondrial function is profiled in a new way by measuring the rates of dye reduction from electrons flowing into and through the electron transport chain from substrates whose oxidation produces NADH or FADH<sub>2</sub>.




**Figure 3.** The MitoPlate S-1 simultaneously assays the metabolic rates of potential NADH or FADH<sub>2</sub> producing substrates.

A1 No Substrate	A2 D-Glucose	A3 Citrovit	A4 D-Glucose-1-PH4	A5 No Substrate	A6 D-Glucose	A7 Citrovit	A8 D-Glucose-1-PH4	A9 No Substrate	A10 D-Glucose	A11 Citrovit	A12 D-Glucose-1-PH4
B1 D-Glucose-1-PH4	B2 D-Glucose-1-PH4	B3 D,L-Glycerol-PH4	B4 L-Lactic Acid	B5 D-Glucose-1-PH4	B6 D-Glucose-1-PH4	B7 D,L-Glycerol-PH4	B8 L-Lactic Acid	B9 D-Glucose-1-PH4	B10 D-Glucose-1-PH4	B11 D,L-Glycerol-PH4	B12 L-Lactic Acid
C1 Pyruvic Acid	C2 Citric Acid	C3 D,L-Succinic Acid	C4 cis-Isocitric Acid	C5 Pyruvic Acid	C6 Citric Acid	C7 D,L-Succinic Acid	C8 cis-Isocitric Acid	C9 Pyruvic Acid	C10 Citric Acid	C11 D,L-Succinic Acid	C12 cis-Isocitric Acid
D1 L-Aspartic Acid	D2 Succinic Acid	D3 Pyruvic Acid	D4 L-Malic Acid	D5 L-Aspartic Acid	D6 Succinic Acid	D7 Pyruvic Acid	D8 L-Malic Acid	D9 L-Aspartic Acid	D10 Succinic Acid	D11 Pyruvic Acid	D12 L-Malic Acid
E1 L-Aspartic Acid	E2 D,L-β-Hydroxybutyric Acid	E3 L-Glutamic Acid	E4 L-Glutamine	E5 L-Aspartic Acid	E6 D,L-β-Hydroxybutyric Acid	E7 L-Glutamic Acid	E8 L-Glutamine	E9 L-Aspartic Acid	E10 D,L-β-Hydroxybutyric Acid	E11 L-Glutamic Acid	E12 L-Glutamine
F1 α-Glc	F2 β-Glc	F3 α-Glc	F4 Fructose	F5 α-Glc	F6 β-Glc	F7 α-Glc	F8 Fructose	F9 α-Glc	F10 β-Glc	F11 α-Glc	F12 Fructose
G1 Malic Acid	G2 Acetyl-L-Carnitine	G3 Carnitine	G4 Palmitoyl-DL-Carnitine	G5 Malic Acid	G6 Acetyl-L-Carnitine	G7 Carnitine	G8 Palmitoyl-DL-Carnitine	G9 Malic Acid	G10 Acetyl-L-Carnitine	G11 Carnitine	G12 Palmitoyl-DL-Carnitine
H1 Pyruvic Acid	H2 L-Aspartic Acid	H3 L-Aspartic Acid	H4 L-Aspartic Acid	H5 Pyruvic Acid	H6 L-Aspartic Acid	H7 L-Aspartic Acid	H8 L-Aspartic Acid	H9 Pyruvic Acid	H10 L-Aspartic Acid	H11 L-Aspartic Acid	H12 L-Aspartic Acid

**Figure 4.** The MitoPlate I-1 simultaneously assays the sensitivity of NADH or FADH<sub>2</sub> producing pathways to 22 mitochondrial inhibitors.



A1 No inhibitor	A2 No inhibitor	A3 No inhibitor	A4 No inhibitor	A5 No inhibitor	A6 No inhibitor	A7 No inhibitor	A8 No inhibitor	A9 No inhibitor	A10 No inhibitor	A11 No inhibitor	A12 No inhibitor
B1 Complex I inhibitor	B2 Complex I inhibitor	B3 Complex I inhibitor	B4 Complex I inhibitor	B5 Complex I inhibitor	B6 Complex I inhibitor	B7 Complex I inhibitor	B8 Complex I inhibitor	B9 Complex I inhibitor	B10 Complex I inhibitor	B11 Complex I inhibitor	B12 Complex I inhibitor
C1 Complex II inhibitor	C2 Complex II inhibitor	C3 Complex II inhibitor	C4 Complex II inhibitor	C5 Complex II inhibitor	C6 Complex II inhibitor	C7 Complex II inhibitor	C8 Complex II inhibitor	C9 Complex II inhibitor	C10 Complex II inhibitor	C11 Complex II inhibitor	C12 Complex II inhibitor
D1 Complex III inhibitor	D2 Complex III inhibitor	D3 Complex III inhibitor	D4 Complex III inhibitor	D5 Complex III inhibitor	D6 Complex III inhibitor	D7 Complex III inhibitor	D8 Complex III inhibitor	D9 Complex III inhibitor	D10 Complex III inhibitor	D11 Complex III inhibitor	D12 Complex III inhibitor
E1 Rotenone	E2 Rotenone	E3 Rotenone	E4 Rotenone	E5 Rotenone	E6 Rotenone	E7 Rotenone	E8 Rotenone	E9 Rotenone	E10 Rotenone	E11 Rotenone	E12 Rotenone
F1 Malonate	F2 Malonate	F3 Malonate	F4 Malonate	F5 Malonate	F6 Malonate	F7 Malonate	F8 Malonate	F9 Malonate	F10 Malonate	F11 Malonate	F12 Malonate
G1 FCCP	G2 FCCP	G3 FCCP	G4 FCCP	G5 FCCP	G6 FCCP	G7 FCCP	G8 FCCP	G9 FCCP	G10 FCCP	G11 FCCP	G12 FCCP
H1 Valproic Acid	H2 Valproic Acid	H3 Valproic Acid	H4 Valproic Acid	H5 Valproic Acid	H6 Valproic Acid	H7 Valproic Acid	H8 Valproic Acid	H9 Valproic Acid	H10 Valproic Acid	H11 Valproic Acid	H12 Valproic Acid
I1 Carbonyl	I2 Carbonyl	I3 Carbonyl	I4 Carbonyl	I5 Carbonyl	I6 Carbonyl	I7 Carbonyl	I8 Carbonyl	I9 Carbonyl	I10 Carbonyl	I11 Carbonyl	I12 Carbonyl
J1 Polystyrene II	J2 Polystyrene II	J3 Polystyrene II	J4 Polystyrene II	J5 Polystyrene II	J6 Polystyrene II	J7 Polystyrene II	J8 Polystyrene II	J9 Polystyrene II	J10 Polystyrene II	J11 Polystyrene II	J12 Polystyrene II

**Figure 5.** Colon and liver cells were assayed for substrate metabolism using the MitoPlate S-1. Four major differences were found in their metabolism.



**Figure 6.** Colon cells were assayed in MitoPlate I-1 with L-malate as the substrate. The cells were sensitive to Complex 1 and 3 inhibitors, but not Complex 2 inhibitors.

**Figure 7.** Colon cells were assayed in MitoPlate I-1 with succinate as the substrate. The cells were sensitive to Complex 2 and 3 inhibitors, but not Complex 1 inhibitors.



## Conclusions

The MitoPlate assay technology enables profiling of mitochondrial function in much greater detail. The MitoPlates have 53 phenotypic assays already dried in the wells, so they can be tested at the same time by simply inoculating with a cell suspension. The assays are colorimetric and can be performed using any kinetic microplate reader. The technology provides a simple and highly sensitive discovery tool for mitochondrial researchers.